
Ovulation Test™

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Resultados / Risultati

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FOR *IN VITRO* USE ONLY ENGLISH

■ CLINICAL APPLICATION:

• DETERMINATION OF THE FERTILE PERIOD

The physical signs of oestrus (vulval swelling, vaginal discharge, or attraction to the male) only provide an approximate indication of ovulation and are not sufficient to determine the optimal timing of mating or insemination. Depending on the bitch, the exact timing of ovulation can be highly variable: 20% of bitches ovulate outside of the period between the 10th and 14th days of heat. Vaginal smears, to confirm whether the bitch is in oestrus or pro-oestrus, combined with progesterone assays, provides the most accurate determination of ovulation to enable scheduling of insemination or mating at the ideal time for optimal results.

The Luteinizing Hormone (LH) spike triggers ovulation, which occurs within 2 days. After ovulation, oocyte maturation takes 2 to 3 days, and oocytes remain viable for 48 to 72 hours. The most fertile period for the bitch is therefore 2 to 5 days after ovulation. Pregnancy will only be possible with mature ovules (2 days after ovulation) and live sperm (5 days on average).

Progesterone level is low just prior to ovulation (between 0 to 1 ng/ml). It begins to increase as LH peak, reaching 6 to 8 ng/ml at the time of ovulation; after ovulation, progesterone level continues to increase over the next few days and remains high during the overall gestation. In the end of gestation, progesterone level decreases, lower than 1 ng/ml. When used in conjunction with all other available informations such as reduced rectal temperature, white milk, ultrasound examination, etc., a progesterone concentration of less than 1 ng/ml is indicative of parturition within 24 hours to 48 hours.

■ PRINCIPLE

Ovulation Test is a rapid, membrane-based, semi-quantitative ELISA

test with visual reading for the determination of progesterone levels in the plasma or serum of bitches. It is essentially used to determine timing of ovulation and therefore to predict the optimal date for insemination or mating. Ovulation Test is an ideal complement to a vaginal smear examination during oestrus monitoring.

■ PROCEDURE

► FOR EACH TEST YOU WILL NEED:

One test cup, one single-use pipette, the bottle of Substrate A, the bottle of Substrate B, the bottle of Wash Solution 1, the bottle of Wash Solution 2, the bottle for mixing the substrates, the bottle of Conjugate, and the colour chart.

Method:

Take the kit out from the refrigerator at least 2 hours before use to allow the reagents to come to room temperature. The kit can be left at room temperature all day long.

Do not mix reagents from different batches.

It is very important to respect the waiting times between each stage of the test (especially steps 1, 3 and 7) to be confident with the accuracy of the results.

1/ With a pipette, add 8 drops of sample (serum or plasma) to the centre of the cup.

Wait 2 minutes.

2/ Add 4 drops of Wash solution 1 to the centre of the cup.

Wait until the liquid has been completely absorbed.

Repeat this step once.

3/ Add 3 drops of conjugate (bottle with red cap) to the centre of the cup.

Wait 1 minute.

4/ Fill the cup to the lower border of the rim with Wash solution 2, i.e. approximately 20 drops.

Wait until the liquid has been completely absorbed.

5/ Prepare a fresh substrate solution using the empty mixing bottle (blue label): add one volume each of substrate A and substrate B (black caps) by filling their droppers to the black line. Homogenize by shaking the freshly prepared substrate. Do not prepare the substrate in advance because it starts to deteriorate after 30 minutes.

6/ Add 4 drops of the freshly prepared substrate (blue label) to the centre of the cup. Discard any unused substrate solution after 30 minutes as it rapidly becomes unstable. Keep the empty vial for subsequent tests.

7/ Wait 9 minutes before reading the result.

Use the colour chart to interpret the result.

■ INTERPRETATION OF RESULTS

The colour corresponds to the progesterone concentration. Use the colour chart supplied with the kit to interpret and record the result.

The following table is a guide to the use of the test for monitoring oestrus prior to insemination in the bitch.

Colour	Progesterone concentration	Interpretation and action*
C1 Bright Blue	0 to 1 ng/ml	Baseline progesterone concentration. Repeat the test every 2 days until a light blue colour is obtained.
C2 Light Blue	1 to 2.5 ng/ml	LH peak. Progesterone level is increasing. Test again every day until the test result becomes C3 (faint blue).
C3 Faint Blue	2.5 to 8 ng/ml	The eggs begin to ripen after ovulation. Mate or inseminate between 1 and 3 days later (5 days with frozen semen).
C4 White	> 8 ng/ml	Mate or Inseminate immediately. Caution: the fertile period might be over.

* The veterinarian should interpret all test results in light of the patient's history, ultrasound examinations, and/or vaginal smears, before making any decisions.

■ RECOMMENDATIONS

• SHELF LIFE:

- 14 months from the manufacturing date written on the box. Store between +2°C and +8°C. Do not freeze.
- The expiry date is printed on each box.

• SAMPLES:

- Collect blood samples into an EDTA or Heparin coated tube or a plain tube. Do not use tubes containing serum separator gel. Do not use whole blood.
- Immediately after collection, invert the tube several times to mix. Centrifuge or allow to clot by leaving at room temperature for 2 hours, then decant the serum or plasma into a clean glass tube. Do not use a haemolysed or lipemic sample (it should be crystal clear).
- For testing within 24 hours: store serum or plasma at room temperature.
- Testing between 24 to 72 hours of sampling: store serum or plasma between +2°C and +8°C (refrigerator).
- Testing more than 72 hours after sampling: store serum or plasma at -20°C

• PRECAUTIONS:

- To ensure correct migration of the sample, place the test device on a horizontal surface.
- To standardize the amount of sample recommended for each test, use only the single-use pipettes supplied in the kit.
- To ensure that the sample and reagent are correctly applied, always hold the pipette or dropper bottle vertically.

• WHEN TO TEST:

- Vaginal cytology can be used to identify the approach of the LH peak and thus when to start testing progesterone levels.
- 1 or 2 tests are usually sufficient if the first day of the onset of

pro-oestrus is known (vaginal bleeding), given that ovulation occurs 2 days after the LH peak (between the 10th and 14th day after the beginning of oestrus).

- The frequency of use of the Ovulation Test depends on each case as a function of the timing of oestrus and the patient's history. In some cases, such as those with a fertility problem or irregular cycles, it is advisable to start testing from the 4th day of pro-oestrus and to repeat the test every 2 days until insemination.

•**PRECAUTIONS:**

- Store the kit in the fridge when not in use.
- Always reseal the plastic bag after removing a test cup.
- **Do not mix test cups or reagents from different kits.**

•**OTHER RECOMMENDATIONS:**

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to assist practitioners in the detection of ovulation and in determining the optimal timing of mating or insemination. All test results should therefore be interpreted in the light of the patient's history, and clinical and gynaecological examination. The definitive diagnosis remains the veterinarian's prerogative and responsibility. Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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